



PDLIM1 Antibody

Product Code	CSB-PA017731GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00151
Immunogen	Human PDLIM1,CLP36
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	PDZ and LIM domain 1;PDLIM1;CLIM1;CLP-36;CLP36;hCLIM1 ;Elfin
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PDLIM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.