



PLA2G12B Antibody

Product Code	CSB-PA018087GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BX93
Immunogen	Human PLA2G12B-Specific
Species Reactivity	Human, Mouse
Tested Applications	ELISA
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	22kd;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLA2G12B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.