



POTED Antibody

Product Code	CSB-PA001758GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YR6
Immunogen	Human POTED-Specific
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	POTE ankyrin domain family, member D; POTED; A26B3; ANKRD21; POTE; POTE-21; POTE21 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POTED
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.