



POU3F1 Antibody

Product Code	CSB-PA018396GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03052
Immunogen	Human POU3F1-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	POU class 3 homeobox 1;POU3F1;OCT6;OTF6;SCIP ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POU3F1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.