



# RMDN3 Antibody

<b>Product Code</b>	CSB-PA008337GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96TC7
<b>Immunogen</b>	Human PTPIP51
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC, IF
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	family with sequence similarity 82, member A2; FAM82A2; FAM82C; FLJ10579; RMD3; hRMD-3; ptpip51 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Apoptosis
<b>Target Names</b>	RMDN3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.