



RMDN3 Antibody

Product Code	CSB-PA008337GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96TC7
Immunogen	Human PTPIP51
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	family with sequence similarity 82, member A2; FAM82A2; FAM82C; FLJ10579; RMD3; hRMD-3; ptpip51 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Apoptosis
Target Names	RMDN3