



DDIT4 Antibody

Product Code	CSB-PA006590GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX09
Immunogen	Human REDD1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	DNA-damage-inducible transcript 4;DDIT4;Dig2;FLJ20500;REDD-1;REDD1;RP11-442H21.1;RTP801 ;
Immunogen Species	Homo sapiens (Human)
Target Names	DDIT4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.