



REXO4 Antibody

Product Code	CSB-PA019587GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZR2
Immunogen	Human REXO4,REX4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	REX4, RNA exonuclease 4 homolog (S. cerevisiae);REX4;REX4;XPMC2;XPMC2H ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	REXO4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.