



RIMS2 Antibody

Product Code	CSB-PA019723GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ26
Immunogen	Human RIMS2-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	regulating synaptic membrane exocytosis 2;RIMS2;DKFZp781A0653;KIAA0751;OBOE;RAB3IP3;RIM2 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RIMS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.