



SENP5 Antibody

Product Code	CSB-PA021011GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HI0
Immunogen	Human SENP5-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SUMO1/sentrin specific peptidase 5;SENP5;DKFZp564O1016;FLJ42398;MGC27076 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SENP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.