



SEPT9 Antibody

Product Code	CSB-PA021033GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHD8
Immunogen	Human SEPT9
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	septin 9;SEPT9;AF17q25;FLJ75490;KIAA0991;MSF;MSF1;NAPB;PNUTL4;SINT1;SeptD1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SEPT9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.