



TGFB2 Antibody

Product Code	CSB-PA023448GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61812
Immunogen	Human TGFB2-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	transforming growth factor, beta 2;TGFB2;MGC116892;TGF-beta2 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	TGFB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.