





NDUFB7 Antibody

Product Code	CSB-PA01594A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17568
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 7 protein (2-137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 7 (Cell adhesion protein SQM1) (Complex I-B18) (CI-B18) (NADH-ubiquinone oxidoreductase B18 subunit), NDUFB7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NDUFB7

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA01594A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA01594A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.