



HLA-DRB1 Antibody

Product Code	CSB-PA17949A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01912
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DRB1-3 chain protein (31-266AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DRB1-3 chain (Clone P2-beta-3) (MHC class II antigen DRB1*3), HLA-DRB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HLA-DRB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.