




polA Antibody

Product Code	CSB-PA16667A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19821
Immunogen	Recombinant Thermus aquaticus DNA polymerase I, thermostable protein (310-832AA)
Raised In	Rabbit
Species Reactivity	Thermus aquaticus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase I, thermostable (EC 2.7.7.7) (Taq polymerase 1), polA, pol1
Immunogen Species	Thermus aquaticus
Research Area	Others
Target Names	polA

Image



Western blot

All lanes: Aquaticus DNA polymerase I, thermostable antibody at 2µg/ml + 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 94 kDa

Observed band size: 94 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.