



polA Antibody, FITC conjugated

Product Code	CSB-PA16667C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19821
Immunogen	Recombinant Thermus aquaticus DNA polymerase I, thermostable protein (310-832AA)
Raised In	Rabbit
Species Reactivity	Thermus aquaticus
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase I, thermostable (EC 2.7.7.7) (Taq polymerase 1), polA, pol1
Immunogen Species	Thermus aquaticus
Research Area	Others
Target Names	polA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.