



# Ifng Antibody

Product Code	CSB-PA15574A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01581
Immunogen	Recombinant Rat Interferon gamma protein (23-156AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon gamma (IFN-gamma), Ifng
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Immunology
Target Names	Ifng