




# TIMP3 Antibody

|                            |                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA11749A0Rb                                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                            |
| <b>Uniprot No.</b>         | P35625                                                                                                   |
| <b>Immunogen</b>           | Recombinant Human Metalloproteinase inhibitor 3 protein (30-208AA)                                       |
| <b>Raised In</b>           | Rabbit                                                                                                   |
| <b>Species Reactivity</b>  | Human                                                                                                    |
| <b>Tested Applications</b> | ELISA, WB; Recommended dilution: WB:1:1000-1:5000                                                        |
| <b>Form</b>                | Liquid                                                                                                   |
| <b>Conjugate</b>           | Non-conjugated                                                                                           |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                         |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                 |
| <b>Isotype</b>             | IgG                                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                               |
| <b>Alias</b>               | Metalloproteinase inhibitor 3 (Protein MIG-5) (Tissue inhibitor of metalloproteinases 3) (TIMP-3), TIMP3 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                     |
| <b>Research Area</b>       | Cell Biology                                                                                             |
| <b>Target Names</b>        | TIMP3                                                                                                    |

## Image



Western blot  
 All lanes: TIMP3 antibody at 2µg/ml  
 Lane 1: EC109 whole cell lysate  
 Lane 2: 293T whole cell lysate  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 25 kDa  
 Observed band size: 25 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.