




ureB Antibody

Product Code	CSB-PA14359A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69996
Immunogen	Recombinant Helicobacter pylori Urease subunit beta protein (1-569AA)
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Helicobacter pylori
Research Area	Signal Transduction
Target Names	ureB

Image



Western Blot

Positive WB detected in: ureB 293T Transfected lysate, 293T non-Transfected lysate

All lanes: ureB antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 90 kDa

Observed band size: 90 kDa