



ureB Antibody, Biotin conjugated

Product Code	CSB-PA14359D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69996
Immunogen	Recombinant Helicobacter pylori Urease subunit beta protein (1-569AA)
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Urease subunit beta (EC 3.5.1.5) (Urea amidohydrolase subunit beta), ureB, hpuB
Immunogen Species	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
Research Area	Signal Transduction
Target Names	ureB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.