



ACTB Antibody, HRP conjugated

Product Code	CSB-PA01629B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Recombinant Human Actin, cytoplasmic 1 protein (102-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin, cytoplasmic 1 (Beta-actin) [Cleaved into: Actin, cytoplasmic 1, N-terminally processed], ACTB
Immunogen Species	Homo sapiens (Human)
Research Area	Isotype/Loading Controls
Target Names	ACTB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.