



LCAT Antibody, FITC conjugated

Product Code	CSB-PA11369C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04180
Immunogen	Recombinant Human Phosphatidylcholine-sterol acyltransferase protein (25-433AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylcholine-sterol acyltransferase (EC 2.3.1.43) (Lecithin-cholesterol acyltransferase) (Phospholipid-cholesterol acyltransferase), LCAT
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LCAT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.