




# IFNG Antibody

<b>Product Code</b>	CSB-PA14709A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8CGS0
<b>Immunogen</b>	Recombinant Guinea Pig Interferon gamma protein (24-166AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Guinea Pig
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	IFN-gamma, IFNG
<b>Immunogen Species</b>	Guinea Pig
<b>Research Area</b>	Others
<b>Target Names</b>	IFNG

## Image



Western blot  
 All lanes: IFNG antibody at 2µg/ml +  
 Recombinant Guinea Pig Interferon gamma  
 protein 10µg  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 20 kDa  
 Observed band size: 20 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.