



AHSG Antibody, FITC conjugated

Product Code	CSB-PA06369C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02765
Immunogen	Recombinant Human Alpha-2-HS glycoprotein protein (20-296AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-2-HS-glycoprotein (Alpha-2-Z-globulin) (Ba-alpha-2-glycoprotein) (Fetuin-A) [Cleaved into: Alpha-2-HS-glycoprotein chain A; Alpha-2-HS-glycoprotein chain B], AHSG, FETUA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	AHSG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.