



CTSF Antibody

Product Code	CSB-PA15279A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBX1
Immunogen	Recombinant Human Cathepsin F protein (273-484AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin F (CATSF) (EC 3.4.22.41), CTSF
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CTSF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.