



Map3k1 Antibody

Product Code	CSB-PA15289A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53349
Immunogen	Recombinant Mouse Mitogen-activated protein kinase kinase kinase 1 protein (1216-1493AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 1 (EC 2.7.11.25) (MAPK/ERK kinase kinase 1) (MEK kinase 1) (MEKK 1), Map3k1, Mekk Mekk1
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Map3k1