



Map3k1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA15289A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P53349 |
| Immunogen | Recombinant Mouse Mitogen-activated protein kinase kinase kinase 1 protein (1216-1493AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Mitogen-activated protein kinase kinase kinase 1 (EC 2.7.11.25) (MAPK/ERK kinase kinase 1) (MEK kinase 1) (MEKK 1), Map3k1, Mekk Mekk1 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Signal Transduction |
| Target Names | Map3k1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |