



Map3k1 Antibody, FITC conjugated

Product Code	CSB-PA15289C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53349
Immunogen	Recombinant Mouse Mitogen-activated protein kinase kinase kinase 1 protein (1216-1493AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 1 (EC 2.7.11.25) (MAPK/ERK kinase kinase 1) (MEK kinase 1) (MEKK 1), Map3k1, Mekk Mekk1
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Map3k1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.