




poxB Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA17967D0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P07003 |
| Immunogen | Recombinant Escherichia coli Pyruvate dehydrogenase [ubiquinone] protein (1-572AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Tested Applications | ELISA, WB; Recommended dilution: WB:1:500-1:5000 |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Pyruvate dehydrogenase [ubiquinone] (EC 1.2.5.1) (Pyruvate oxidase) (POX) [Cleaved into: Alpha-peptide], poxB |
| Immunogen Species | Escherichia coli |
| Research Area | Others |
| Target Names | poxB |

Image



Western Blot

Positive WB detected in: Rosseta bacteria lysate at 40µg, 20µg, 10µg, 5µg

All lanes: poxB antibody at 0.25µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 63 kDa

Observed band size: 63 kDa