




poxB Antibody, Biotin conjugated

Product Code	CSB-PA17967D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07003
Immunogen	Recombinant Escherichia coli Pyruvate dehydrogenase [ubiquinone] protein (1-572AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyruvate dehydrogenase [ubiquinone] (EC 1.2.5.1) (Pyruvate oxidase) (POX) [Cleaved into: Alpha-peptide], poxB
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	poxB

Image



Western Blot

Positive WB detected in: Rosseta bacteria lysate at 40µg, 20µg, 10µg, 5µg
 All lanes: poxB antibody at 0.25µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 63 kDa
 Observed band size: 63 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.