



IL18BP Antibody, Biotin conjugated

Product Code	CSB-PA16739D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95998
Immunogen	Recombinant Human Interleukin-18-binding protein (31-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-18-binding protein (IL-18BP) (Tadekinig-alfa), IL18BP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL18BP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.