



clpS Antibody, Biotin conjugated

Product Code	CSB-PA16649D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A8Q6
Immunogen	Recombinant Escherichia coli ATP-dependent Clp protease adapter protein ClpS protein (1-105AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent Clp protease adapter protein ClpS, clpS, yljA
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	clpS