



Nos3 Antibody, HRP conjugated

Product Code	CSB-PA16579B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q62600
Immunogen	Recombinant Rat Nitric oxide synthase, endothelial protein (519-690AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nitric oxide synthase, endothelial (EC 1.14.13.39) (Constitutive NOS) (cNOS) (EC-NOS) (Endothelial NOS) (eNOS) (NOS type III) (NOSIII), Nos3
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Nos3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.