



Hp Antibody, Biotin conjugated

Product Code	CSB-PA17069D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06866
Immunogen	Recombinant Rat Haptoglobin protein (19-346AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Haptoglobin (Liver regeneration-related protein LRRG173) (Zonulin) [Cleaved into: Haptoglobin alpha chain; Haptoglobin beta chain], Hp
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Hp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.