



mpt64 Antibody, HRP conjugated

Product Code	CSB-PA14949B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WIN9
Immunogen	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunogenic protein MPT64 (Antigen MPT64), mpt64
Immunogen Species	Mycobacterium tuberculosis
Research Area	Others
Target Names	mpt64
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.