




PLAU Antibody

Product Code	CSB-PA14319A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00749
Immunogen	Recombinant Human Urokinase-type plasminogen activator protein (21-173AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Urokinase-type plasminogen activator (U-plasminogen activator) (uPA) (EC 3.4.21.73) [Cleaved into: Urokinase-type plasminogen activator long chain A; Urokinase-type plasminogen activator short chain A; Urokinase-type plasminogen activator chain B], PLAU
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLAU

Image



Western blot

All lanes: PLAU antibody at 1.8µg/ml

Lane 1: Mouse heart tissue


Lane 2: Mouse skeletal muscle tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 49, 47 kDa

Observed band size: 49 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA14319A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.