




CDH12 Antibody

Product Code	CSB-PA13109A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55289
Immunogen	Recombinant Human Cadherin-12 protein (56-605AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin-12 (Brain cadherin) (BR-cadherin) (Neural type cadherin 2) (N-cadherin 2), CDH12
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CDH12

Image



Western blot

All lanes: CDH12 antibody at 2µg/ml

Lane 1: EC109 whole cell lysate


Lane 2: 293T whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 89, 85 kDa

Observed band size: 85 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA13109A0Rb at dilution of 1:50



Immunofluorescent analysis of PC-3 cells cells using CSB-PA13109A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.