




Vegfa Antibody

Product Code	CSB-PA08249A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16612
Immunogen	Recombinant Rat Vascular endothelial growth factor A protein (27-232AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vascular endothelial growth factor A (VEGF-A) (Vascular permeability factor) (VPF), Vegfa, Vegf
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Vegfa

Image



Western Blot

Positive WB detected in: Rat liver tissue


All lanes: Vegfa antibody at 4.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 26, 23, 21, 18 kDa

Observed band size: 26 kDa



Immunohistochemistry of paraffin-embedded mouse lung tissue using CSB-PA08249A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.