




TNC Antibody

Product Code	CSB-PA11504A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24821
Immunogen	Recombinant Human Tenascin protein (1888-2201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tenascin (TN) (Cytotactin) (GMEM) (GP 150-225) (Glioma-associated-extracellular matrix antigen) (Hexabrachion) (JI) (Myotendinous antigen) (Neuronection) (Tenascin-C) (TN-C), TNC, HXB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TNC

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA11504A0Rb at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.