





ABCB1 Antibody

Product Code	CSB-PA11737A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08183
Immunogen	Recombinant Human Multidrug resistance protein 1 protein (236-297AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multidrug resistance protein 1 (EC 3.6.3.44) (ATP-binding cassette sub-family B member 1) (P-glycoprotein 1) (CD antigen CD243), ABCB1, MDR1 PGY1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABCB1


Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA11737A0Rb at dilution of 1:50



Immunofluorescent analysis of HepG2 cells using CSB-PA11737A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA11737A0Rb at dilution of 1:100