




MAP1S Antibody

Product Code	CSB-PA13549A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q66K74
Immunogen	Recombinant Human Microtubule-associated protein 1S protein (794-1053AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated protein 1S (MAP-1S) (BPY2-interacting protein 1) (Microtubule-associated protein 8) (Variable charge Y chromosome 2-interacting protein 1) (VCY2-interacting protein 1) (VCY2IP-1) [Cleaved into: MAP1S heavy chain; MAP1S light chain], MAP1S, BPY2IP1 C19orf5 MAP8 VCY2IP1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAP1S

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA13549A0Rb at diluton of 1:50

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.