



SERPINE1 Antibody, HRP conjugated

Product Code	CSB-PA11239B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05121
Immunogen	Recombinant Human Plasminogen activator inhibitor 1 protein (24-398AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasminogen activator inhibitor 1 (PAI) (PAI-1) (Endothelial plasminogen activator inhibitor) (Serpine E1), SERPINE1, PAI1 PLANH1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPINE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.