



SERPINE1 Antibody, HRP conjugated

Product Code	CSB-PA11239B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05121
Immunogen	Recombinant Human Plasminogen activator inhibitor 1 protein (24-398AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasminogen activator inhibitor 1 (PAI) (PAI-1) (Endothelial plasminogen activator inhibitor) (Serpine E1), SERPINE1, PAI1 PLANH1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPINE1