



AMBP Antibody, FITC conjugated

Product Code	CSB-PA10009C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02760
Immunogen	Recombinant Human Protein AMBP protein (20-203AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein AMBP [Cleaved into: Alpha-1-microglobulin (Protein HC) (Alpha-1 microglycoprotein) (Complex-forming glycoprotein heterogeneous in charge); Inter-alpha-trypsin inhibitor light chain (ITI-LC) (Bikunin) (EDC1) (HI-30) (Uronic-acid-rich protein); Trypstatin], AMBP, HCP ITIL
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	AMBP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.