



Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 6 (CEACAM6)

Product Code	CSB-MP005166HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40199
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	KLTIES TPFNVAEGKE VLLLAHNLPQ NRIGYSWYKG ERVDGNSLIV GYVIGTQQAT PGPAYSGRET IYPNASLLIQ NVTQNDTGFY TLQVIKSDLV NEEATGQFHV YPELPKPSIS SNNSNPVEDK DAVAFTCEPE VQNTTYLWWV NGQSLPVS PR LQLSNGNMTL TLLSVKRND A GSYECEIQNP ASANRSDPVT LNVLYGPDGP TISPSKANYR PGENLNLSCH AASNPPAQYS WFINGTFQQS TQELFIPNIT VNNSGSYMCQ AHNSATGLNR TTVTMITVSG
Source	Mammalian cell
Target Names	CEACAM6
Protein Names	Recommended name: Carcinoembryonic antigen-related cell adhesion molecule 6 Alternative name(s): Non-specific crossreacting antigen Normal cross-reacting antigen CD_antigen= CD66c
Expression Region	35-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.