



INHBB Antibody, HRP conjugated

Product Code	CSB-PA13929B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09529
Immunogen	Recombinant Human Inhibin beta B chain protein (295-405AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibin beta B chain (Activin beta-B chain), INHBB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INHBB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.