



KLK10 Antibody, Biotin conjugated

Product Code	CSB-PA14039D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43240
Immunogen	Recombinant Human Kallikrein-10 protein (31-274AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-10 (EC 3.4.21.-) (Normal epithelial cell-specific 1) (Protease serine-like 1), KLK10, NES1 PRSSL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KLK10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.