



ORM1 Antibody, Biotin conjugated

Product Code	CSB-PA05709D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02763
Immunogen	Recombinant Human Alpha-1-acid glycoprotein 1 protein (19-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-1-acid glycoprotein 1 (AGP 1) (Orosomuroid-1) (OMD 1), ORM1, AGP1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ORM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.