



IdcC Antibody, FITC conjugated

Product Code	CSB-PA16349C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52095
Immunogen	Recombinant Escherichia coli Constitutive lysine decarboxylase protein (1-713AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Constitutive lysine decarboxylase (LDC) (EC 4.1.1.18), IdcC, Idc IdcH
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	IdcC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.