




HPRT1 Antibody

Product Code	CSB-PA15827A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00492
Immunogen	Recombinant Human Hypoxanthine-guanine phosphoribosyltransferase protein (2-218AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hypoxanthine-guanine phosphoribosyltransferase (HGPRT) (HGPRTase) (EC 2.4.2.8), HPRT1, HPRT
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HPRT1

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA15827A0Rb at dilution of 1:20

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.